


coloris™

Identification Clips for PIPETMAN®

EN



www.gilson.com/contactus

Gilson, Inc.

3000 Parmenter Street • PO Box 620027
Middleton, WI 53562 USA
T 608-836-1551 or 800-445-7661 • F 608-831-4451

Gilson S.A.S.

19, avenue des Entrepreneurs BP 145 • F-95400
Villiers-le-Bel, France
T +33 (0) 1 34 29 50 00 • F +33 (0) 1 34 29 50 20

LT801278/I © 06/2017 Gilson Inc. All rights reserved. All trademarks are the property of Gilson Inc. and its subsidiaries.

 **GILSON®**

 **GILSON®**

Coloris™ identification clips are compatible with:

PIPETMAN Classic™ P2 to P1000,
PIPETMAN® Fixed F2 to F1000,
PIPETMAN Neo® P2N to P1000N,
PIPETMAN® G P2G to P1000G,
PIPETMAN® L P2L to P10mLL,
PIPETMAN® L Fixed F1L to F5000L.

ORDERING REFERENCE	DESCRIPTION
F161301	10 units, mixed colours
F161302	10 units, red
F161303	10 units, yellow
F161304	10 units, green
F161305	10 units, blue
F161306	10 units, white

Personalize your PIPETMAN with a touch of color!



- Avoids accidental exchanges.
- Available in five colors.
- Clear and clean.

Write your name or the lab's name on the clip with a permanent ink pen. Wipe the clip with isopropanol to remove the ink.

Made of POM (Polyacetal)